

Cartridges Designed for a Diverse Range of Needs



SYSTEM SPECIFICATIONS

Energy Type	HIFU (High Intensity Focused Ultrasound)
Power	0.1 - 3.0J (0.1J / Step)
Rating	100-240V~, 50/60Hz
Weight	10kg
Dimension	450 (L) X 550 (W) X 270 (H)mm (Including Holder)
Cartridge	B2-6.0 (2MHz / 6.0mm) B2-9.0 (2MHz / 9.0mm) B2-13.0 (2MHz / 13.0mm)



UTB02BB2018



CONFIDENTLY FIT with

Ulfit

Macro Focused Circular Ultrasound

Body Shaping
Experience that
Satisfies Worldwide

317,000+

Ulfit

Body Slimming & Tightening

Recognized across the globe for its advanced technology of MFCU, the ULFIT achieves faster and more effective treatments designed exclusively for all body profiles and areas targeted. Hence, a proven device that maximizes patient comfort and satisfaction in between every treatment.



A Technology with a
Leading Edge


MFCU Technology

The Industry's First Macro Focused Circular Ultrasound

The circular design of the multi-depth body cartridges is equipped with MFCU technology that delivers wider coverage as well as equal and precise distribution of high-density ultrasound energy into the focal region.




Transducer Technology



Larger Focal Width

Wider transducer yields less pain and precise distribution upon making contact.



Optional Pain Control

Users may select between "G" (General) and "P" (Pain control) modes to dramatically reduce pain levels.

80%

Reduced Pain

The ULFIT is engineered to deliver treatments with 80% less pain.

Perfect Treatment for the Body

Circular Motion Treatment

The ULFIT's body applicator distributes ultrasound energy into large focal regions both evenly and in quick procession.



Convenient Treatment

- One-Touch Continuous Procedure
- Easy Grip & Lightweight



Application




- Axilla
- Arms
- Bra Line
- Abdomen
- Flanks
- Love Handles
- Below Buttocks
- Outer Thighs
- Inner Thighs
- Knees

Accelerated Shot Speed


The high-peak power capabilities of the body cartridge achieve faster treatment times for equally impressive body outcomes.

88:00 / 100 SHOTS




Superior Results

BODY TREATMENT



80% of 15 patients achieved a significant reduction of abdominal fat after 1 treatment with the ULFIT.

Degree of Improvement



Target Patients: 15 individuals in the 20-50 age range with concerns about abdominal circumference. **Testing Period:** 6 weeks. **Testing Methods:** ULFIT 6.0 / 9.0mm Body Cartridges **No. of Treatments:** 1. Both baseline and post treatment measurements of each patient at the 2-week mark were evaluated by three clinicians. **Assessment:** Fair: 1-9mm / Good: 10-24mm / Excellent: 25mm and over.

PATIENT SATISFACTION


89% of 100 participants have reported significant reductions to their body fat as well as impressive changes after receiving treatments with the ULFIT.

Patient Ratio of Improvement



Target Patients: 100 individuals with concerns about body fat in multiple areas. **Testing Methods:** ULFIT 6.0 / 9.0mm Body Cartridges **No. of Treatments:** 1-6. Both baseline and post treatment photos of the patients were evaluated every 1 week post treatment.


Age 38 / F



Age 35 / F




Age 35 / F



Age 34 / F



Age 45 / F



Age 28 / F

